Gebrauchsanweisung Operating Instructions Instrucciones de uso Manuel d'instructions Istruzioni per l'uso

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REUSABLE ROTATING TOOLS

- Contents:**
- Safety precautions
 - Proper use
 - Indication
 - Handling
 - Sterilisation
 - Storage
 - Basic information
 - Working aids

Safety precautions

MEDICON tools are manufactured of high-quality materials. They are developed, designed and produced with the greatest care. Constant further development and strictest quality monitoring are standard procedure at MEDICON.

Only proper use of these quality tools enables the best possible working results and a long service life of the tools. To ensure this, the following operating instructions and safety precautions must always be observed and complied with.

Improper use and use for other than the intended purpose may result in damage to bones and tissue, to premature wear or to destruction of the tools and danger to the user, the patient or third parties.

It is the responsibility of the user to check the product on his/her own responsibility for suitability and the possibility of use for the intended purpose prior to its use.

These operating instructions accompany the following rotating tools:

- Rose-head burr in all shapes
- Hobbing cutter
- Lindemann burr
- Neuro burr
- Craniotomy burr
- Cement burr
- Diamond drill and grinding wheel
- Twist drill
- Spinal burr
- Perforation drill
- Circular saw blades incl. carrier
- Dental drills and cutters

No trephines of any kind are covered by the descriptions in these operating instructions.

Proper use

Rotating tools must be used only in accordance with these operating instructions and the operating instructions for the motor systems.

Before use it must be ensured that the tools employed are in faultless technical and sterile condition. Generally, tools which are solid or damaged must not be used.

Use tools according to the specifications of the motor manufacturer and check that they are securely clamped.

Use tools in accordance with their intended purpose. Comply with the speed specifications for the individual tools and motor systems.

Indication

MEDICON tools have been developed in accordance with their intended purpose for processing bone, skin, soft tissue and cement. Tools with various designs are available for specific purposes.

The surgeon is responsible for selecting the appropriate tool.

Speeds

Excessive speeds must be avoided. They can lead to premature wear and can result in increased heat development and the danger of thermal necroses.

In extreme cases the tool can even break. The maximum permissible speed must always be observed for safety reasons.

HANDLING

Contact force

It is essential for excessive tool contact forces to be avoided. Excessive contact forces can lead to damage to the tool in the form of chipping of the cutting edges and bending. At the same time this can result in increased heat development and the danger of thermal necroses.

In addition, when using grinding wheels, an excessive contact force can cause diamond grains to break away, which can in turn result in increased heat development.

In extreme cases the tool may break.

It is essential to ensure that the tool is not tilted at an angle during surgery. Never use the tool as a lever.

Cooling/Rinsing

When using rotating tools, sufficient cooling or rinsing must always be provided. Inadequate cooling or rinsing will lead to immediate clogging of the cutting edges, and therefore to increased heat development. (Danger of thermal necroses) In addition, the service life of the instruments may be considerably shortened.

Dull and damaged tools

Check the tools before use with a magnifying glass (6-10 x magnification) for bluntness and damage, in which the following should be checked:

- Chipping of the cutting edges
- Dislodging of diamond grains
- Bare metal areas
- Damage to the shank
- Bent tools

REUSABLE ROTATING TOOLS

Blunt and damaged instruments must never be used. Failure to observe these precautions may result in tool breakage and danger to the patient and the user!

Cleaning/disinfection

The rotating tools must be cleaned, disinfected and sterilised before initial use and after each additional use.

Tools must be laid in a suitable disinfectant solution immediately following the operation, and blood, secretion, tissue and bone residues must be removed immediately. (**danger of drying on**)

Equipment

Cleaning/disinfection device (RDG) (Miele with VarioTD program) or similar, adjustable to the parameters of the Vario TD program (See diagram in fig. 1).

Alkaline cleaner (please observe manufacturer's specifications)

1. Thoroughly rinse off tools with running water immediately prior to machine-processing to prevent cleaning agent and disinfectant residues from getting into the machine.
2. Place the tools in a suitable tool carrier or screen basket.
3. Position the tool carrier in the RDG so that the rinsing spray strikes the tool directly.
4. Start the Vario-TD program (see diagram in fig. 1) incl. thermal disinfection. The thermal disinfection must be carried out according to the AO value and national regulations (EN/ISO 15883).

The process operates according to a fixed pattern:

- 4 minutes pre-wash with cold water

- Draining
- 6 minutes pre-wash with alkaline cleaner
- Draining
- 3 minutes neutralisation with warm tap water (> 40°C)
- Draining
- 2 minutes intermediate rinse with warm tap water (> 40°C)
- Draining

5. Removal of residual moisture from tools with filtered compressed air which does not lead to recontamination.

If visible soiling is still present on the tool following machine processing, cleaning and disinfection must be repeated.

Sterilization

Only cleaned and disinfected tools may be sterilized.

Only the sterilization methods listed may be used for sterilization.

- Fractionated vacuum method with sufficient product drying
- Processing must be carried out according to a validated method.

Steam sterilization in the fractionated vacuum method in a sterilizer in accordance with DIN EN 13060 or DIN EN 285, validated according to DIN EN ISO 17665-1 or DIN EN ISO 14937.

**Sterilization temperature = 134°C
Minimum dwell time = 5 min.**

A temperature of over 180°C must always be avoided, as this will lead to reduced tool hardness and shorten the service life.

To avoid spotting and corrosion, the steam must be free of additives. The recommended limits for additives for feed water and condensed steam are specified in DIN EN 285.

When sterilizing several tools, the maximum loading of the sterilizer must not be exceeded. (See manufacturer's specifications)

Storage

The processed rotating tools must be stored in germ-tight packaging, dry and protected from dust in suitable and regularly disinfected cabinets or similar storage locations.

Basic information

Comply with the applicable regulations for processing medical products in the country concerned.

The manufacturer must ensure that the processing methods listed above are suitable for processing the tools, which are not yet in use but have been removed from the sterile packaging. The operating authority is responsible for ensuring that the processing actually carried out in the reprocessing device achieve the desired results with the equipment, materials and personnel employed. Validation and routine monitoring of the procedure are required for this purpose. Any deviation from the process described here (e.g. use of other process chemicals) should be carefully evaluated by the processor for their effectiveness and possible negative consequences.

MEDICON eG, as the company marketing these products, does not assume any liability for direct or consequential damage resulting from improper use or handling, in particular as a result of failure to comply with the operating instructions or as a result of incorrect care or maintenance.

Working aids	Tool after use	Preparation/Time	Visual Inspection/Cleaning	Preparation	Procedure Temperature/Time	Provision
	Disinfection		Cleaning/Care		Sterilization	Storage
	Immediately place in instrument disinfection and cleaning solution	Use DGHM-recognised disinfectant for rotating instruments ①	Check instruments for possible defects and residues with a magnifying glass (10x magnification)	Place in instrument carrier, Steri-Container (DIN 58946) ②	Steam (autoclave) 134°C / 3,5 - 5 min 121°C / 15 - 20 min ③	In Steri-Container (DIN 58946) ④ ⑤
	and/or cleaning in ultrasonic bath	Time: According to manufacturers specifications	Remove residues with brush, rinse off and blow off with compressed air	Seal individual instruments in sterile material packaging (DIN 58953)		or sterile material packaging (DIN 58953) ④ ⑤
	Important: Thermal disinfection is not suitable for instrument disinfection and cleaning.	Use DGHM-listed method Time: According to manufacturers specifications	Blow off with compressed air, rinse off and blow off again			

① Observe information on materials when selecting the disinfection method, the disinfection preparation and the sterilization method.

② The sterilization process must be guaranteed when the instrument carrier and the sterile material packages are selected.

③ Minimum soaking times (DIN 58946, Part 1). The operating times may vary depending on the manufacturer.

④ The type of sterile material package must be based on the sterilisation method used (DIN 58953).

⑤ Guide values for storing sterile materials (DIN 58953).

Important: Adhering disinfectant can behave aggressively under pressure and temperature and can result in damage to the instruments.

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