

MASTER SERIES

A Micro Bond System

Metal Brackets

Maxillary	Torq	Ang	Rot	M-D	R/L	.018		.022		
						No hook	With hook	No hook	With hook	
Central	+12	+5	0	Diagonal micro with .020 x .020 VS	.050	R L	390-0080 390-0081		390-0090 390-0091	
Lateral	+9	+8	0	Diagonal micro with .020 x .020 VS	.050	R L	390-0082 390-0083		390-0092 390-0093	
Cuspid	-2	+12	0	Diagonal micro with .020 x .020 VS	.050	R L	390-0084 390-0085		390-0094 390-0095	
1st and 2nd Bicuspid	-10	0	0	Diagonal micro, hook gingival with .020 x .020 VS	.050	R/L	390-0089	390-0089B	390-0099	390-0099B

Mandibular										
Anteriors	0	0	0	Beveled diagonal micro with .020 x .020 VS	.050	R/L	390-0086		390-0096	
Cuspid	-11	+8	0	Beveled diagonal micro with .020 x .020 VS	.050	R L	390-0087 390-0088		390-0097 390-0098	
1st and 2nd Bicuspid	-15	0	0	Diagonal micro, hook gingival with .020 x .020 VS	.050	R/L	390-0189	390-0189B	390-0199	390-0199B


BRACKET ID
Color coded by quadrant



HGW: high gingival wing
D-G: distal/gingival
M-G: mesial/gingival
VS: vertical slot

PERMANENT MARKS

Maxillary Central, lateral, cuspid - dimple gingival
1st and 2nd bicuspid - ID line occlusal
Mandibular 1st and 2nd bicuspid - ID line occlusal



MASTER SERIES

A Friction Free System

Maxillary	Torq	Ang	Rot	M-D	R/L	.018		.022		
						No hook	With hook	No hook	With hook	
Central	+12	+5	0	Single	.120	R L	590-0001 590-0002		590-0021 590-0022	
Lateral	+8	+9	0	Single	.120	R L	590-0003 590-0004		590-0023 590-0024	
Cuspid	-2	+10	0	Single, hook D-G	.120	R L	590-0005 590-0006	590-0005B 590-0006B	590-0025 590-0026	590-0025B 590-0026B
1st and 2nd Bicuspid	-7	0	0	Single, no HGW	.120	R/L	590-0007		590-0027	

Mandibular										
Anteriors	0	0	0	Single, bevel occlusal	.105	R/L	590-0008		590-0028	
Cuspid	-11	+7	0	Single, bevel occlusal, hook D-G	.120	R L	590-0009 590-0010	590-0009B 590-0010B	590-0029 590-0030	590-0029B 590-0030B
1st Bicuspid	-17	0	0	Single, bevel occlusal	.120	R/L	590-0011		590-0031	
2nd Bicuspid	-22	0	0	Single, bevel occlusal, no HGW	.120	R/L	590-0012		590-0032	

BRACKET ID
Color coded by quadrant



HGW: high gingival wing
D-G: distal/gingival
M-G: mesial/gingival
VS: vertical slot

PERMANENT MARKS

Maxillary Central, lateral, cuspid - dimple gingival
1st and 2nd bicuspid - dimple gingival
Mandibular Cuspid - dimple gingival
2nd bicuspid - dimple gingival

