

Mini Tweed

Maxillary	Torq	Ang	Rot	M-D	R/L	.018		.022		
						No hook	With hook	No hook	With hook	
Central	0	0	0	.100	R/L	380-0201	For hook option add 000-013 to part number	380-0221	For hook option add 000-013 to part number	
Lateral	0	0	0	.100	R/L	380-0202		380-0222		
Cuspid	0	0	0	.100	R/L	380-0203		380-0223		
Bicuspid	0	0	0	.080	R/L	380-0214		380-0234		
Mandibular										
Anterior	0	0	0	.100	R/L	380-0208	For hook option add 000-013 to part number	380-0228	For hook option add 000-013 to part number	
Anterior	0	0	0	.080	R/L	380-0218		380-0238		
Cuspid	0	0	0	.100	R/L	380-0203		380-0223		
Bicuspid	0	0	0	.080	R/L	380-0214		380-0234		

Recommended wire diameters: .017 x .022, .018 x .025, .019 x .025, .0215 x .028

BRACKET ID
NOT color coded

