

# MINI MASTER SERIES

## A Rinchuse\* Bio-Dimensional Vertical Slot System

Maxillary	Torq	Ang	Rot	M-D	R/L	.018		.022	
						No hook	With hook	No hook	With hook
Central	+17	+4	0	.140	R L	390-6101 390-6102			
Lateral	+10	+8	0	Hook D-G	R L	390-6103 390-6104	390-6103B 390-6104B		
Cuspid	0	+8	0	Hook D-G	R L			390-6125 390-6126	390-6125B 390-6126B
1st and 2nd Bicuspid	-7	0	0	Hook D-G	R L				390-6128BD 390-6129BD
OPTION	-7	0	0	OFFSET DB bicuspid tube with L hook	2 mm R L				098-1406C 098-1407C

### Mandibular

Central	-6	0	0	Bevel occlusal	.100	R/L	390-6118		
Lateral	-6	0	0	Bevel occlusal, hook D-G	.100	R/L		390-6118B	
Cuspid	-7	+3	0	Bevel occlusal, hook D-G	.122	R L		390-6159 390-6160	390-6159B 390-6160B
1st Bicuspid	-12	+2	0	Bevel occlusal, hook D-G	.120	R L		390-6131 390-6132	390-6131B 390-6132B
2nd Bicuspid	-17	+2	0	Bevel occlusal, hook D-G	.120	R L		390-6133 390-6134	390-6133B 390-6134B
OPTION	-17	0	0	OFFSET DB bicuspid tube with L hook	2 mm	R L			098-1414C 098-1415C

**BRACKET ID**  
Color coded by quadrant




**HGW:** high gingival wing  
**D-G:** distal/gingival  
**M-G:** mesial/gingival  
**VS:** vertical slot

**OFFSET BICUSPID OPTIONS:** Gingivally offset bicuspid brackets offer a larger bonding base covering more area on the occlusal one third of the buccal surface and increasing bond strength.

#### PERMANENT MARKS

**Maxillary** Central, lateral, cuspid - dimple D-G  
1st and 2nd bicuspid - ID line occlusal

**Mandibular** Cuspid - dimple D-G  
2nd bicuspid - ID lines occlusal



\*Dr. Don Rinchuse and Dr. Dan Rinchuse are proud to endorse Rinchuse prescription brackets from American Orthodontics.