



# SURGICAL AND RECTANGULAR ARCH HOOK


- Prevent failure of bracket attachments to teeth
- 2 mm cast with ball or "L" hook
- Split hooks can be placed without removing arch wire
- Provide an efficient means of ligation
- Use sliding hook crimping plier for strongest crimp
- Surgical hooks may be used for intermaxillary fixation following surgery




## SURGICAL HOOK


	Fluted Slide on	.018 x .025	851-620CF
		.021 x .025	851-621CF
			10/pkg


	2 mm Slide on	.018 x .025	851-620C
		.021 x .025	851-621C
			10/pkg


	2 mm Split	.018 x .025	851-622C
		.021 x .025	851-623C
			10/pkg

## ARCH HOOK

	Fluted Slide on	.018 x .025	851-632CF
		.021 x .025	851-633CF
			10/pkg

	2 mm Slide on	.018 x .025	851-632C
		.021 x .025	851-633C
			10/pkg

	2 mm Split Right	.018 x .025	851-630CR
		.021 x .025	851-631CR
			10/pkg

	2 mm Split Left	.018 x .025	851-630CL
		.021 x .025	851-631CL
			10/pkg



## Sliding Hook Crimping Plier

- Two grooved inserts maintain a firm grip on the tube
- Used with surgical hooks or rectangular arch hooks

100-107 For surgical hooks (shown)  
100-107A For arch hooks

## CRIMPABLE

Round Crimpable	.020	
2 mm Arch Hook	.022	851-701
	.025	851-702
	.030	851-703
	.036	851-704
	.040	851-705
	.045	851-706
	.051	851-707
Wall thickness: .006		10/pkg



Round Crimpable	.020	851-740X
1 mm Stop	.022	851-741X
	.025	851-742X
	.030	851-743X
	.036	851-744X
	.040	851-745X
	.045	851-746X
	.051	851-747X
Wall thickness: .006		100/pkg



Round Crimpable	.020	851-740
2 mm Stop	.022	851-741
	.025	851-742
	.030	851-743
	.036	851-744
	.040	851-745
	.045	851-746
	.051	851-747
Wall thickness: .006		100 /pkg



Round Split Crimpable	.030	851-634
2 mm Stop		
Wall thickness: .015		50/pkg



Cast Rectangular	.018 x .025	851-635C
Split Crimpable	.021 x .028	851-636C
2 mm Stop		
Wall thickness: .015		10/pkg



Rectangular	.018 x .025	851-612
Split Crimpable	.021 x .028	851-613
2 mm Stop		
Wall thickness: .015		100/pkg



## Crimpable

- Ideal for light wires
- Use dull hard wire cutter for placement
- Can place split crimpable without removing arch wire

## Weldable

- Welded to arch wire prior to ligation

## WELDABLE

Round Weldable	.022	851-720
2 mm Arch Hook	.036	851-723
	.040	851-724
	.045	851-725
	.051	851-726
Wall thickness: .015		10/pkg




Round Weldable	.036	851-763
2 mm Stop	.040	851-764
	.045	851-765
	.051	851-766
Wall thickness: .015		100/pkg



Rectangular Weldable	.018 x .025	851-614
1.25 mm Stop	.021 x .028	851-615
Wall thickness: .015		100/pkg



Rectangular Weldable	.018 x .025	851-610CF
2 mm Fluted Stop	.021 x .028	851-611CF
		100/pkg
2 mm Stop	.018 x .025	851-610
	.021 x .028	851-611
Wall thickness: .015		100/pkg



Arch hook

Round arch stop

Rectangular arch stop

## MONKEY HOOK



292-186 Bondable coil eyelet

### 851-640

- .018 linkable S-shaped hook
- Directs impacted teeth into position
- Assists in the correction of severely rotated teeth
- May be hung from a bondable coil eyelet #292-186
- Can be used in multiples to create a chain
- 10 pieces per package

## ARCH LOCK

- Removable lock secures to arch wires from .012 (0.30 mm) up to .022 x .028 (0.56 x 0.71 mm) and retainer wire up to .030 (0.76 mm)
- Pre-assembled unit allows placement in one piece


### **KIT** 851-650 Arch Lock Kit Without Hooks

- 10 Arch locks without hooks  
2 Hex keys



### **KIT** 851-651 Arch Lock Kit With Hooks

- 10 Arch locks with hooks  
2 Hex keys



### Reorder Components

- 851-652 10 locks without hooks  
851-653 10 locks with hooks  
851-655 5 hex keys