


EAGLEWEAR

Safety Release Headgear System

Safety Release Spring Modules reduce dangerous "slingshot" rebound injuries while providing measurable, consistent force levels.

Colorful **EAGLEWEAR** products help keep patients enthusiastic about wearing their headgear.



852-812

A EAGLEWEAR High Pull Headgear Strap

- Dual loop design allows for built in size adjustment
- Teal on black
- Packaged individually

852-828

B EAGLEWEAR Slimline High Pull Release Module

- For use with Eaglewear high-pull headgear strap only
- 10 modules (5 pair) per package

C EAGLEWEAR Cervical Neck Strap

- Dual loop design for small or large neck size patients
- Washable material for a clean looking strap
- Use with Safety Release Spring Module
- 10 straps per package

852-745 Pink on black

852-746 Purple on black

852-747 Teal on black

852-748 Assorted

EAGLEWEAR Pouches

- Zip closure and patient I.D. label
- 10 pouches per package

852-741 Pink

852-742 Purple

852-743 Teal

852-744 Assorted

D Safety Release Spring Module

- Use with safety cervical neck straps or high-pull headgear straps
- Reduces risk of "slingshot" rebound of the facebow
- Available in three force levels, with color-coded attachment straps
- Modules available in dark navy blue and white
- Ten modules (5 pair) per package

Spring force

852-820 Light - 16 oz. - Navy/clear strap

852-821 Medium - 20 oz. - Navy/yellow strap

852-825 Medium - 20 oz. - White/white strap

852-822 Heavy - 23 oz. - Navy/blue strap

852-815

Safety Cervical Neck Strap

- Dark navy blue
- Smooth seams for patient comfort
- 10 straps per package

852-811

Safety High Pull Headgear Strap

- Dark navy blue to blend with spring release modules
- 1 strap per package

J-Hook Headcap Assembly

- Headcap and J-Hook pre-assembled for your convenience
- Color coded cheek pads for quick size selection
- J-Hooks easily adaptable for fit and comfort
- Use Posted arches with J-Hook assembly
- Use #000-133 Extra Heavy elastics with J-Hook assembly

852-730 Small - beige on black

852-731 Medium - blue on black

852-732 Large - black on black

Use posted arches and #000-133
extra heavy elastics with the
J-Hook assembly

**Ortho Band Neck Pads and Straps**

852-750 7" blue satin backed foam rubber neck pad - pkg of 10

852-752 5-1/2" blue satin backed foam rubber neck pad - pkg of 10

852-758 1-1/2" wide blue satin backed foam rubber neck pad - 5 yard roll

852-756 3/4" blue adjustable elastic strap - pkg of 10





852-754 Safety chain - pkg of 1



NOTE: Elastic straps will not be sold without a safety chain.



300 SERIES FACE BOWS



WARNING: Use only with safety cervical headgear strap or safety high pull headgear strap.

- .072 spring tempered outer bow
- Inner bow reinforced with inconel tubing
- Outer bow designed to avoid frontal cheek area
- Centerline permanently stamped: .045=single line .051=double line

 Standard style	Standard Outer Bow		Short Outer Bow		Long Outer Bow		
	.045 inner	852-311	852-315	852-318	.051 inner	852-312	852-316
 Standard Outer Bow style							

 Loop style	Size 2 90 mm		Size 3 97 mm		Size 4 104 mm		Size 5 111 mm		Size 6 118 mm	
	.045 inner	852-302	852-303	852-304	852-305	.051 inner	852-322	852-323	852-324	852-325
 Loop style Sizes marked: .045 = black .051 = red										



 Loop style with short outer bow	Size 2 90 mm		Size 3 97 mm		Size 4 104 mm		Size 5 111 mm		Size 6 118 mm	
	.045 inner	852-372	852-373	852-374	852-375	852-376				
 Loop style with short outer bow Sizes marked: .045 = black										

 Adjustable loop with molar offset bend	Size 2 90 mm		Size 3 97 mm		Size 4 104 mm		Size 5 111 mm		
	.045 inner	852-345	852-347	852-349	852-351	.051 inner	852-346	852-348	852-350
 Adjustable loop with molar offset bend Sizes marked: .045 = black .051 = red									

Size is determined by measuring along the circumference of the inner bow from distal of the offset bend to distal of the offset bend.

440 SERIES FACE BOWS

- .062 outer bow
- Inconel sleeve extended on outer bow

 Root style	Size 1 83 mm		Size 2 90 mm		Size 3 97 mm		Size 4 104 mm		Size 5 111 mm	
	.045 inner	852-431	852-432	852-433	852-434	.051 inner	852-430	852-436	852-437	852-438
 Root style Sizes marked: .045 = black .051 = red										